



# Video Capture Library for LabVIEW

December 2015

## Version 3.4

---

**HAMAMATSU**

Homepage Address <http://www.hamamatsu.com>

HAMAMATSU PHOTONICS K.K., Systems Division

812 Joko-cho, Hamamatsu City, 431-3196, Japan, Telephone: (81)53-431-0124, Fax: (81)53-435-1574, E-mail: [export@sys.hpk.co.jp](mailto:export@sys.hpk.co.jp)

U.S.A. and Canada: Hamamatsu Photonic Systems: 360 Foothill Road, Bridgewater, N.J. 08807-0910, U.S.A., Telephone: (1)908-231-1116, Fax: (1)908-231-0852, E-mail: [usa@hamamatsu.com](mailto:usa@hamamatsu.com)

Germany: Hamamatsu Photonics Deutschland GmbH: Arzbergerstr. 10, D-82211 Herrsching am Ammersee, Germany, Telephone: (49)8152-375-0, Fax: (49)8152-2658, E-mail: [info@hamamatsu.de](mailto:info@hamamatsu.de)

France: Hamamatsu Photonics France S.A.R.L.: 8, Rue du Saule Trapu, Parc du Moulin de Massy, 91882 Massy Cedex, France, Telephone: (33)1 69 53 71 00, Fax: (33)1 69 53 71 10, E-mail: [infos@hamamatsu.fr](mailto:infos@hamamatsu.fr)

United Kingdom: Hamamatsu Photonics UK Limited: 2 Howard Court, 10 Tewin Road Welwyn Garden City Hertfordshire AL7 1BW U.K., Telephone: (44)0 1707-294888, Fax: (44)0 1701-325777, E-mail: [info@hamamatsu.co.uk](mailto:info@hamamatsu.co.uk)

North Europe: Hamamatsu Photonics Norden AB: Smidesvägen 12, SE-171-41 Solna, Sweden, Telephone: (46)8-509-031-00, Fax: (46)8-509-031-01, E-mail: [info@hamamatsu.se](mailto:info@hamamatsu.se)

Italy: Hamamatsu Photonics Italia S.R.L.: Strada della Moia, 1/E 20020 Arese (Milano), Italy, Telephone: (39)02-935 81 733, Fax: (39)02-935 81 741, E-mail: [info@hamamatsu.it](mailto:info@hamamatsu.it)

<b>Overview .....</b>	<b>4</b>
<b>System Requirements.....</b>	<b>4</b>
<b>What's New In Version 3 .....</b>	<b>4</b>
<b>Installation .....</b>	<b>5</b>
Hamamatsu Video Capture Installation.....	5
Camera Driver Installation .....	5
<b>Standard Functions.....</b>	<b>6</b>
TM_INITIALIZE .....	6
TM_DEINITIALIZE .....	6
TM_OPENCAMERA .....	6
TM_CLOSECAMERA .....	7
TM_GETCAMERAINFO .....	7
TM_SETPARAMETER.....	7
TM_GETPARAMETER .....	7
TM_GETRANGE .....	8
TM_GETFUNCTIONLIST .....	8
TM_SETAREA .....	9
TM_GETAREA.....	9
TM_SETINPUTTRIGGER.....	10
TM_SETOUTPUTTRIGGER .....	10
TM_FIRETRIGGER.....	11
TM_PREPARECAPTURE .....	11
TM_UNPREPARECAPTURE .....	11
TM_STARTCAPTURE.....	11
TM_STOPCAPTURE .....	12
TM_WAITNEXTFRAME .....	12
TM_GETCAPTUREINFO .....	12
TM_GETFRAME.....	12
<b>DSP Functions.....</b>	<b>13</b>
TM_RECURSIVEFILTER .....	13
TM_FRAMEAVERAGING.....	13
TM_SETBACKGROUNDFRAME .....	13
TM_SETSHADINGFRAME .....	13
TM_STOREBACKGROUNDFRAME .....	14
TM_STORESHADINGFRAME.....	14
<b>Hard Disk Recording Functions.....</b>	<b>15</b>
TM_STARTRECORDER .....	15
TM_STOPRECORDER.....	15
TM_GETRECORDERSTATUS .....	15
TM_OPENDCIMGFILE.....	16
TM_CLOSEDCIMGFILE.....	16
TM_GETDCIMGFRAME .....	16
<b>Samples .....</b>	<b>19</b>
Single Frame Acquisition .....	19
Preview Capture .....	19
Adjusting Camera Settings.....	20

Function Enumeration .....21

Recording to a DCIMG File.....22

Reading a DCIMG File .....23

Converting Intensity values to Electron values.....23

Setting Region of Interest .....24

Getting Region of Interest .....24

# Overview

The Hamamatsu Video Capture Library software is a collection of VIs designed to control and gather image data from a Hamamatsu camera from within LabVIEW. Because this is built using Hamamatsu's DCAMAPI, applications created with this library will be able to detect and use many of the Hamamatsu line of cameras. And because it is made for LabVIEW, a programmer will be able to create large virtual instruments in a relatively short amount of time.

This library can be used with IMAQ Vision from National Instruments. IMAQ Vision is a powerful image processing library that is made for LabVIEW. IMAQ Vision is not required to use this software.

It is understood that the user is familiar with the Hamamatsu camera(s) and its functions. If you are unfamiliar with some of the capabilities of the camera, please refer to your camera manual. You will find technical specifications and installation instructions. If you encounter a problem with your camera that you cannot solve, please contact our support group.

It is also understood that the user is familiar with National Instruments LabVIEW. If you need further assistants with LabVIEW, please refer to your LabVIEW manual or contact National Instruments.

## System Requirements

- Microsoft Windows 7 or higher
- Hamamatsu camera with compatible DCAMAPI module
- National Instruments LabVIEW 2011 or higher

## What's New In Version 3


This version of the Video Capture Library has updated the VIs to include additional functions that provide you additional controls. Preparing a capture session and starting a capture session are now two different functions. Waiting for a new frame to arrive and collecting the data into LabVIEW are also two different functions. These changes have been made to allow for maximum performance so that you can collect all of your frame data from some of our more demanding camera systems.

Also included are functions for the new Hard Disk Recording feature. These functions allow you to capture data and write the frames directly to the hard disk into DCIMG format. This method is ideal for capturing a large number of frames in the most efficient way possible. These DCIMG files can be read back later with the functions in this library.

# Installation

## Hamamatsu Video Capture Installation

1. Be sure that LabVIEW is not running. If so, close all LabVIEW windows. Also be sure that your camera drivers have been properly installed. If not, please see the Camera Driver Installation below.
2. Run the SETUP.EXE program that comes with the installation software. For the 64-bit version of LabVIEW, run the SETUP\_X64.EXE program instead to install the 64-bit version.
3. Follow all of the prompts during the installation process.
4. If you have LabVIEW 2011 or newer, the setup will install all of the proper files.
5. To access your new video capture functions, just follow the following diagram to locate the library.




## Camera Driver Installation

1. If you do not have a DCAM-API driver disk, you can download the camera driver from <http://www.dcamapi.com> and unzip the contents to a new directory.
2. If the CD does not automatically run, launch SETUP.EXE from the root directory.
3. Select and click the button from the menu of the interface that your camera is using.
4. Follow the on-screen instructions to install the driver you selected.


# Standard Functions

## TM\_INITIALIZE




This is the first function that you need to run before you can use any camera control functions in the library. Its purpose is to initialize all of the required resources for the video capture library and find all supported cameras on the system. Any camera you wish to use with your system must be connected to your system and powered on before this function is called. When this function returns, it will report the number of supported cameras that were found, the version of this capture library and the DCAM version, and it will also provide the model number, serial number, and firmware version of each camera found.

## TM\_DEINITIALIZE




If all connected devices have been closed, this will cleanly free all resources used to by the video capture library. However, if there are still some cameras opened, this will not force close those devices. It is recommended that you close all devices before calling this function.

## TM\_OPENCAMERA




This will activate a camera for operation at the specified Camera Index and will return a handle for that camera as well as some details about the camera. These details include the model number, serial number, and firmware version of the camera. This handle must be used by the other camera control functions to determine the camera to be controlled. The index values start at 0, so if there are 3 cameras available on the system, the valid input for Camera Index would be 0, 1, and 2. This function can only be used to activate cameras that were found by TM\_INITIALIZE. If the camera is already open in another VI or was not closed properly, this function will fail. If you set Force Open to TRUE, the function will attempt to close any open session with the camera before opening the camera.

## TM\_CLOSECAMERA




This will close the camera associated with the handle and free resources that it used. Once this function is called, the handle will become invalid and cannot be used anymore. The camera can still be opened again with TM\_OPENCAMERA.

## TM\_GETCAMERAINFO




This function will return information about the camera. The possible items to get information are the camera model, camera ID, camera bus, vendor, driver version, camera version, module version, and DCAMAPI version. This function can also be used on unopened cameras by using the camera index instead of the camera handle.

## TM\_SETPARAMETER




This function can adjust the camera settings. These settings may include gain, offset, exposure time, binning, etc. Not all cameras are capable of every listed function. While you can adjust some settings during a capture, settings that can change the image resolution should not be changed during a capture. For cameras with multiple views, View will select which will be set. Please refer to your camera manual to determine the capabilities of your camera.

## TM\_GETPARAMETER




This function gets the current value of a selected camera setting. Not all cameras are capable of every function. For cameras with multiple views, View will select which value to read.

## TM\_GETRANGE



This function returns the range of the input parameter. The values you will receive will include minimum, maximum, stepping, and default values. If the parameter is a mode function, it will have associated text values and the Text Values output will return true. If that value returns true, then Value Names will contain an array of text values for this function along with their associated number values.


## TM\_GETFUNCTIONLIST



The Hamamatsu camera line is expanding with an ever growing number of functions. This VI will output an array which contains a list of the supported functions available to the camera. This array can be connected directly to the StringsAndValues[] property of a menu ring. The list of supported functions will vary depending on the camera. Additional information of these functions can be obtained from TM\_GETRANGE. And the function values can be controlled with TM\_GETPARAMETER and TM\_SETPARAMETER.




## TM\_SETAREA




Some Hamamatsu cameras are capable of selecting a specific sub-region of the camera to capture image data from. With this function, you will be able to select the area you wish you capture through the Sub Region cluster input. Because of the limitations of the different cameras, the values you enter may not be the same values that were set. It is important to use TM\_GETAREA to see the values that were set to the camera. This should not be called while capturing images.

## TM\_GETAREA



This function will return the value of the current Sub Region. This function will also return the current minimum, maximum, and increment values of each variable. These limits can change whenever you set a new sub-array or binning. It is always a good idea to call this function to get the current valid limits.

## TM\_SETINPUTTRIGGER




This function will allow you to set how the camera will trigger for new images. Parameters contains four settings you can adjust. Mode will allow you to set one of the following trigger modes

- **Internal** - This will set the camera to do its own triggering. This will also produce the fastest frame rate.
- **External edge** – The camera will wait for an external trigger signal to go active.
- **External level** – The camera will wait for an external trigger signal to go active and continue exposing until that signal switches to inactive.
- **Software** – The camera will wait for TM\_FIRETRIGGER to be called.
- **Fast Repetition** – This is a special mode. Not available for all cameras.
- **TDI** – This is a special mode. Not available for all cameras.
- **Trigger Start** – The camera will wait for an external trigger signal to go active and then proceeds to work in internal trigger mode.
- **Sync Readout** – The camera will wait for the first trigger to start an exposure. All of the following triggers will stop the current exposure and start the next one.
- **Master Pulse** – This is a special mode. Not available for all cameras.

Polarity will allow you to set the external trigger to negative or positive. Times is used with Sync Readout. Delay determines the length of time until the camera will react to a trigger. This is only available on some cameras.

## TM\_SETOUTPUTTRIGGER




This function will allow you to set the programmable output triggers of your camera.

Parameters contains four settings that you can adjust. Mode will set the type of output.

- **Low** – Sets the output to low regardless of what the camera is doing.
- **Global Exposure** – Sets the output active when all of the pixels are exposing at the same time.
- **Programmable** – Sets the output to active for a specified period after a specified delay.
- **Trigger Ready** – Sets the output to active when the camera is ready to receive another input trigger.
- **High** – Sets the output to high regardless of what the camera is doing.


Polarity determines which direction is active. Period sets the length of time the output will remain active during Programmable trigger mode. Delay sets the length of time the output will wait before setting the output to active during Programmable trigger mode. The index option determines which output trigger these settings will be applied to.

## TM\_FIRETRIGGER




If your camera is currently in software trigger mode, no images will be acquired until a software trigger is fired. When a software trigger is fired with this function, the camera will begin exposing a new image. This function is very useful for synchronizing the camera with application events.

## TM\_PREPARECAPTURE




This function prepares the driver and camera for capture. There are two inputs available to configure your capture function. CAPTURE MODE allows you to adjust the behavior of the capturing cycle. In Snap mode, the camera will capture the specified number of frames and will stop at the end of the cycle. In Sequence mode, the camera will collect data continuously without stopping at the end of the cycle. Once it reaches the end of a cycle of frames, it will continue on the first frame. FRAMES allow you to adjust the number of frames of images you wish to capture. In Sequence mode, setting frames is unnecessary and will automatically be set to 3. This function will also return an Image Type value which is used by IMAQ Vision to determine the image data type. This function does not start the camera to capture images.

## TM\_UNPREPARECAPTURE




This function will free all resources created with TM\_PREPARECAPTURE. This step is important if you wish to change certain camera settings such as the binning or region of interest. If the camera is capturing data, you must stop it first with TM\_STOPCAPTURE before calling this function.

## TM\_STARTCAPTURE




This begins capturing of image data from the camera to your computer. You will need to have called TM\_PREPARECAPTURE to setup the driver for capture before calling this function.

## TM\_STOPCAPTURE




This function will stop the current capture session. The current capture mode and internal buffers create will remain intact so you will be able to start a new capture session immediately with TM\_STARTCAPTURE.

## TM\_WAITNEXTFRAME




This function will wait for the next frame to become available. When a frame arrives, this function will return and you are able to work with the data. By default, the timeout is set to 5000 ms. The timeout value can be adjusted by setting the timeout input with a new value. If you set a timeout of -1, this function will wait indefinitely until a new frame arrives. The frame index is the buffer index location where the data was stored. Frame count is the number of frames captured since calling TM\_STARTCAPTURE.

## TM\_GETCAPTUREINFO



This function will return the current frame index and frame count of the capture sequence. The frame index is the buffer index location where the data is stored. Frame count is the number of frames captured since calling TM\_STARTCAPTURE.

## TM\_GETFRAME




This function will retrieve the frame in the buffer at the specified frame index. This function requires an array input for the data buffer created by TM\_PREPARECAPTURE. If the output data is 16 bit, the Data 16 output pin will be valid and contain the data in a U16 array. If the output data is 8 bit, the Data 8 output pin will be valid and it will contain the data in a U8 array. Some of the newer cameras provide electron conversion information. If the current camera provides this information, it will be included in the Electron Info output. This conversion info is specific to this frame of data and may be inaccurate if applied to a different frame of data. This is the formula to convert the intensity value to electrons.

$$e^- = (intensity - offset) * coefficient$$

## DSP Functions


Some newer Hamamatsu cameras have imaging processing functions built into the camera hardware. Since these tools have been placed in the hardware, you will not need to use computer resources to perform some of the more common image processing tasks. Please check your camera manual to determine if your camera has these capabilities.

### TM\_RECURSIVEFILTER




This DSP function allows you enable recursive filtering. Each output image of the camera will be the average of the number of frames specified by the Frames input. You can disable this function by setting Frames to 1.

### TM\_FRAMEAVERAGING




This DSP function allows you enable frame averaging. Each output image of the camera will be the average of the number of frames specified by the Frames input. You can disable this function by setting Frames to 1.

### TM\_SETBACKGROUNDFRAME




This DSP function allows you to enable a background image to be subtracted from the output data. You must specify the background frame you wish to subtract with the Frame input.

### TM\_SETSHADINGFRAME




This DSP function allows you to enable a shading correction image to be applied to the output data. You must specify the correction frame you wish to apply with the Frame input.

## TM\_STOREBACKGROUNDFRAME



This DSP function will store a frame of data into the camera memory frame specified by the input Frame for background subtraction. This function will fail if the camera is capturing images.

## TM\_STORESHADINGFRAME




This DSP function will store a frame of data into the camera memory frame specified by the input Frame for shading correction. This function will fail if the camera is capturing images.

# Hard Disk Recording Functions


New to the Hamamatsu Video Capture Library software are the hard disk recorder functions. These functions allow you to capture image data directly to disk. This is especially useful if you are capturing a large amount of data but do not have enough RAM to store the data. The output data is written to a DCIMG file which is a proprietary Hamamatsu image stack format. In order to view images from a DCIMG file, you can use the DCIMG reader functions available with this capture library or you can use any other program designed to read DCIMG files such as HClmage. Please keep in mind that the storage disk is inherently slower than RAM. We recommend using 4 SSDs linked via a RAID0 configuration. This will provide a drive fast enough to handle the data output of our more demanding cameras.

## TM\_STARTRECORDER




This function will enable the hard disk recorder. The recorder will write captured images directly to disk to the specified DCIMG file. This recorder will write up to the specified number of frames. Once that maximum number of frames have been reached, the recording will stop. This function must be called after TM\_PREPARECAPTURE.

## TM\_STOPRECORDER



This function will end the current recorder session and close the DCIMG file. Once the DCIMG file is closed, it cannot be reopened for recording of more images. However, it may be opened by any DCIMG file reader to read all of the data that has been stored.

## TM\_GETRECORDERSTATUS




This function will provide the current index of the newest frame recorded as well as the total number of frames recorded.

## TM\_OPENDCIMGFILE




This function will open a DCIMG file for reading. It will return a file handle which will be used for all other DCIMG related functions. This file handle is different from the camera handle and should not be used for camera control functions. It will also return the total number of frames that are available in the DCIMG file. Once the DCIMG file is opened with this function, the data can be read by TM\_GETDCIMGFRAME.

## TM\_CLOSEDCIMGFILE



This function will release a DCIMG file and free any resources used for it. This will make the file handle invalid and unusable by the other functions.


## TM\_GETDCIMGFRAME



This function will retrieve the frame data from a DCIMG file specified by the file handle and the frame index. The file handle can be retrieved by TM\_OPENDCIMGFILE. The Time Sec output is the timestamp in seconds of the selected frame relative to when the first frame in the file was captured.




## TM\_WRITEMETADATA



This function will allow you to write metadata to a DCIMG file currently open for recording. The Handle is the generated from TM\_STARTRECORDER. With Type, you can specify if this is metadata for a specific frame or for the entire DCIMG file. Both types can be written to the DCIMG file. If you set Type to Frame, then Index will specify which frame the metadata will be written to. Be aware that once TM\_STOPRECORDER is called, the DCIMG handle is destroyed and you will no longer be able to write metadata to the file.

## TM\_READMETADATA



This function will let you read the metadata that is stored in the DCIMG file specified by the file handle generated from TM\_OPENDCIMGFILE. With Type, you can specify if this is metadata for a specific frame or for the entire DCIMG file. Both types can be written to the DCIMG file. If you set Type to Frame, then Index will specify which frame the metadata will be written to.

## Error Code

When using this library, you may encounter an error when running a LabVIEW program. If the error was generated from one of the functions supplied in this library, then you will receive the error code as well as the particular VI that caused the error. This information will provide a clue into determining why the program failed and how it could be corrected. The following is a list of possible error codes that you may encounter.

**TMERR\_ABORT** – The operation has been aborted.

**TMERR\_BUSY** – The camera is currently busy and cannot execute your request.

**TMERR\_FAILEDOPENRECFILE** – Failed to open the specified DCIMG file.

**TMERR\_FAILEDREADDATA** – Failed to read from the DCIMG file.

**TMERR\_FAILEDWRITEDATA** – Failed to write to the DCIMG file.

**TMERR\_FAILOPENCAMERA** – The system was unable to active the camera.

**TMERR\_FAILREADCAMERA** – The camera could not read the command.

**TMERR\_FAILWRITECAMERA** – The camera could not change the parameter.

**TMERR\_INVALIDCAMERA** – The selected device is invalid.

**TMERR\_INVALIDHANDLE** – The specified camera handle is invalid.

**TMERR\_INVALIDPARAM** – The specified parameter is invalid.

**TMERR\_INVALIDPROPERTYID** – The property ID is invalid.

**TMERR\_INVALIDRECHANDLE** – The DCIMG file handle is invalid.

**TMERR\_INVALIDSUBARRAY** – The subarray setting is invalid.

**TMERR\_INVALIDVALUE** – The specified value is invalid.

**TMERR\_LOSTFRAME** – A frame has been lost during transfer from the camera.

**TMERR\_NOCAMERA** – No camera is available.

**TMERR\_NOMEMORY** – Your computer does not have enough memory for your request.

**TMERR\_NORESOURCE** – Your computer does not have enough resources for your request.

**TMERR\_NOTREADY** – The camera is not ready for the command.

**TMERR\_NOTWRITABLE** – The property does not have a writable value.

**TMERR\_NOTREADABLE** – The property does not have a readable value.

**TMERR\_NOTSUPPORT** – The function is not supported by the camera.


**TMERR\_OUTOFRANGE** – The value is out of range.

**TMERR\_TIMEOUT** – The driver has timed out while waiting for the frame data.

**TMERR\_UNKNOWNERROR** – An unknown error has occurred.

# Samples

## Single Frame Acquisition




In this example, we are initializing the library by calling `TM_INITIALIZE` to find all cameras in the system. If it finds at least one camera, we call `TM_OPENCAMERA` with a camera index of 0 to activate the first camera that is found. If successful the function will return a handle to the camera. This handle will be used for all of the other functions.

To acquire a single image, `TM_PREPARECAPTURE` needs to be called to setup the capture. We set the capture mode to Snap and we set only 1 frame. Then we call `TM_STARTCAPTURE` to begin the capture session. `TM_WAITNEXTFRAME` is called to wait for a new frame to be available. Once that function returns, we call `TM_GETFRAME`. We take the Frame Index output from `TM_WAITFRAME` to determine which frame we are going to retrieve from the capture library. The output 2D array will contain the data for the image. We then take this data and send it to the display.


To clean up the acquisition, we first call `TM_STOPCAPTURE` to stop the capture process. Even though this is a snap of 1 image, this is still a necessary step. Then we call `TM_UNPREPARECAPTURE` to free the internal image buffers. Then `TM_CLOSECAMERA` is called to release the camera and all of the resources that were allocated for it. And finally `TM_DEINITIALIZE` is called to unload the capture library.

## Preview Capture



Setting up for a preview is very similar to acquiring a single image. After opening the camera with `TM_OPENCAMERA`, we need to start a sequence mode capture with `TM_STARTCAPTURE`. For inputs, you will only need to set capture mode to Sequence as the frame buffer will automatically be set for you. After the capture is started, we can enter our while loop so we can continuously capture images. `TM_WAITNEXTFRAME` will wait for the next frame. `TM_GETFRAME` will retrieve the next frame from the buffer.


## Adjusting Camera Settings




In this example, we use TM\_OPENCAMERA to create a handle to the Hamamatsu camera. We use TM\_SETPARAMETER to adjust binning mode. For the Parameter input, we use BINNING. And for the Value In input, we use 4 to specify that we want 4x4 binning. We use TM\_SETPARAMETER again, but this time to adjust the gain. For the Parameter input, we use CONTRAST GAIN. For the Value In input, we use 100 to specify that we want to have a gain of 100. Now that we have the settings that we desire, we can then prepare a capture session. Some functions can be adjusted during a capture, however functions that can change the image resolution cannot be changed during a capture. Please refer to your camera manual to determine the capabilities of your camera.

## Function Enumeration


If you are creating a VI capable of handling different models of Hamamatsu cameras, you may choose to enumerate the available functions of the cameras. If you use TM\_GETFUNCTIONLIST with TM\_GETRANGE you will be able to see what functions are available to the camera.



In the above example, the output of TM\_GETFUNCTIONLIST is being used to set the StringsAndValues[] property of the menu ring Function List. This will populate the menu ring with all of the functions available to this camera and associate them with the proper ID value. You may also use the dynamically created Function List menu ring with TM\_GETRANGE. In addition, if the function is a mode function, you can get the available list of mode values for the function using TM\_GETRANGE. Just as in TM\_GETNFUNCTIONLIST, this output can be used to set the StringsAndValues[] property of a menu ring.




The dynamically created Function List and/or Value List menu rings can be used with TM\_SETPARAMETER to control the camera. Be sure that you are using the correct Value List for the function that you are trying to control.




## Recording to a DCIMG File

It is possible to acquire data from the camera and record it directly to the hard drive. Some Hamamatsu cameras are very demanding on the computer that you may not have enough system resources to acquire all of your frames through the standard method. Recording directly to the hard drive removes the need for the data to be processed in LabVIEW first thus allowing more resources to be used to collect the data.



To start a recording, you must first prepare a capture session with TM\_PREPARECAPTURE. You can either snap for the number of frames you wish to record, or set a sequence. Depending on the number of frames you are recording it may be best to set a sequence since this will only create 3 image buffers. If you are recording a large number of frames, you may not have enough RAM. After that function is called, you may now call TM\_STARTRECORDER with the file path and the number of frames to record. Now when images are being captured, they will write directly to the file that you specified.




While it is capturing data, in order to determine the number of frames that were recorded, you will need to use TM\_GETRECORDERSTATUS. This will provide the total number of frames recorded. This will help you to determine if the recording is completed.



Once the recording is done, you should call TM\_STOPRECORDER. This will close the DCIMG file and allows it to be opened by a DCIMG reader. This function can be called anytime during capture.

## Reading a DCIMG File

If you have data that you wish to read from a DCIMG file, you can use the functions available with this capture library. These functions are separate from the main camera control functions and can be used without a camera connected to the computer.



To read a DCIMG file, you will first need to open the file with `TM_OPENDCIMGFILE` and provide the complete file name and path. This will give you the file handle as well as the number of frames available in the DCIMG file. The file handle must be used with the other reader functions.


Once the file is open, you can begin reading the images that are contained in the file by using `TM_GETDCIMGFRAME`. This function requires the file handle provided by `TM_OPENDCIMGFILE` and a frame index. The output will be a 2D array of the image data.

When you are finished with the file, you can close the file with `TM_CLOSEDCIMGFILE`. This will free any resources used to open the file as well as allow any other DCIMG reader to open that file.

## Converting Intensity values to Electron values

Some of the newer camera are capable of providing electron conversion numbers through `TM_GETFRAME`. You can use these numbers to convert the intensity value to electron. Here is the formula to convert the intensity value to electrons


$$e^- = (intensity - offset) * coefficient$$



In this example, the entire image is converted from intensity values to electron values. The data will be floating point (SGL) values and may have negative numbers.


## Setting Region of Interest

The region of interest can be setup in different ways. To set the region, create a cluster of four 32-bit numbers and use this cluster for the Sub Region terminal of TM\_SETAREA. These values are horizontal offset, vertical offset, horizontal size, and vertical size respectfully. You can set the region of interest when the capture session has not been prepared yet.



## Getting Region of Interest

To get the current region of interest or to determine the acceptable values for TM\_SETAREA, you can use TM\_GETAREA. Each of the four outputs are a cluster of four 32-bit values that represent the limits of that variable. In this example, the different values are being used to set the limits of the sub-region controls.





# Parameter Definitions

Some of the capabilities of the Hamamatsu cameras are mode functions. These modes can be changed with this software, however the names of these modes are not always defined. Below is a list of the different mode functions and the definition of their values. Not all cameras are capable of each function or mode. You can get a listing of the camera functions and their value options with the TM\_GETFUNCTIONLIST and TM\_GETRANGE functions. Please refer to your camera manual for more information.

## Sensor Mode

Area	1
Line	3
TDI	4
Progressive (Lightsheet)	12
Split View (W-View)	14

## CCD Mode

Normal CCD	1
EM CCD	2

## Binning

1x1	1
2x2	2
4x4	4
8x8	8
16x16	16

## Light Mode

Low Light Mode	1
High Light Mode	2

## High Dynamic Range Mode

Off	1
On	2

## Sensor Cooler Fan

Off	1
On	2

## Mechanical Shutter

Auto	1
Close	2
Open	3

## Trigger Connector

Interface	1
BNC	2
Multi	3

## Output Trigger Source

Exposure	1
Readout End	2

VSync	3
HSync	4
Trigger	6

### **Trigger Global Exposure**

None	1
Always	2
Delayed	3
Emulate	4
Global Reset	5

### **First Trigger Behavior**

Start Exposure	1
Start Readout	2

### **Readout Direction**

Forward	1
Backward	2
By Trigger	3